

Where are the Johnny Inventors of this world?

A good idea that goes to waste in a drawer due to lack of time or money is a help to nobody.

The inventor's market is looking for good ideas, product improvements and constructive criticism. Of course, we are constantly producing more innovative products and improving our traditional ones.

"But what about all the good ideas of our users that lie dormant?" we have asked ourselves. Now we are asking you directly: You have probably cursed at one time or another, because a tool was good in your daily work but not perfect. Help us to make it better.

With the inventor's market, we want to promote your good ideas and change them into cash for you and us. If you are interested and have one or more ideas for improvements, call us at.

Products developed from customer ideas




Novartos



close-it **MAXI**



Pharma Spoon



ProfiSampler



Flow Indicators



UniSampler with telescopic suction tube

for sampling of liquids from precisely predetermined depths. Telescopic suction tube can be varied at will to a maximum length of 100 cm

Operation

1. Set the desired sampling depth by moving the suction lance.
2. Attach the sampling bottle to the adapter.
3. Lower the suction lance into the medium.
4. Create a vacuum using the handpump and collect the sample.
5. Remove the filled sampling bottle, close and label it.

- PP suction tube, PVC suction hose
- Sample bottles available in 4 sizes

UniSampler with telescopic suction tube

	Price/Unit of (X) units	Article No.	
incl. 5 sample bottles 250 ml and 500 ml each, hand-driven vacuum pump and adapter for sample bottles, 20 seal-it bottle seals	(1) 147,00	5612-1	N4

	Sample bottle 100 ml	(50)	0,78	0319-100	N3
	Sample bottle 250 ml	(25)	1,05	0319-250	N3
	Sample bottle 500 ml	(25)	1,74	0319-500	N3
	Sample bottle 1000 ml	(25)	2,87	0319-1000	N3
	Transport case	(10)	67,00	5330-1	N4