



MicroSampler

useful for taking small samples, sterilizable



1. Pushing in lance **2. Opening head** **3. Close head, pull back lance** **4. Emptying MicroSampler**

The **MicroSampler** is specially designed for removal of the smallest quantities of powders.

Sample quantities of 0.2 ml, 0.5 ml, 1.0 ml, 2.0 ml, 5.0 ml and 10 ml can be taken at varying depths from sacks or containers in a clean, fast and simple way.

The head of the MicroSampler can be exchanged. This enables heads with different volumes to be used.

- Exact sampling point can be chosen
- Different volumes and lengths available
- Conforms to GMP/FDA
- Head removable, easy cleaning
- Made of stainless steel 316
- Special volumes available

Operation

1. Insert MicroSampler into powder.
2. Open the volume head at the target spot by pushing the buttons; the sample flows in.
3. Close the volume head; the sample is collected.
4. Withdraw the lance, open the volume head again, and remove the sample.

Please order lance and head separately! 3 years autoclave 121°

MicroSampler, 12.5 mm dia., without head

Length cm	Tube-diameter mm	Price/Unit	Article No.	
55	12,5	156,00	5307-1055	N4
85	12,5	173,00	5307-1085	N4
120	12,5	188,00	5307-10120	N4

MicroSampler, 25 mm dia., without head

55	25	203,00	5307-2555	N4
85	25	219,00	5307-2585	N4
120	25	234,00	5307-25120	N4

autoclave 121°

Heads fit MicroSampler with 12.5 mm dia.

Tube-diameter mm	Volume ml	Price/Unit	Article No.	
12,5	0,2	46,70	5307-2	N4
12,5	0,5	46,70	5307-5	N4
12,5	1,0	46,70	5307-10	N4

Heads fit MicroSampler with 25 mm dia.

25	2,0	63,00	5307-20	N4
25	5,0	63,00	5307-50	N4
25	10,0	63,00	5307-100	N4



Do you need heads with special volumes?

We'd be glad to help you.

Please enquire!

close-it Quality Control
Contrôle Qualité
Qualitäts-Kontrolle

Datum/Date Name/Nom

close-it

Control seal ensures the tight sealing of punctured sacks/ Bags.